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#### 2. Semi-dry western blot transfer

This protocol describes the procedure for transferring protein from polyacrylamide gel onto membrane using semi-dry transfer method. (Note: It is recommended to transfer very large proteins by wet transfer other than semi-dry transfer ).

0. Transfer buffer composition : 25 mM Tris, 192 mM Glycine, 20% methanol, pH 8,3 (3,03 g Tris, 14,4 g glycine, ddH<sub>2</sub>O, 200 mL methanol or ethanol replacement, qsp 1L ddH<sub>2</sub>O). Optional : add SDS 0,04% final.

1. Cut 12 sheets of whatman 3MM the same size as the PAGE gel. Soak whatman filter paper in transfer buffer for several minutes to 10 min

2.

<p>Nitrocellulose membranes should be handled carefully and with gloves. Nitrocellulose membranes can be brittle and may crack or tear, especially when dry. <u>Cut to the size of the SDS-PAGE gel.</u> Slowly lower the membrane into 20% (v/v) ethanol and agitate briefly. The membrane will become translucent as it wets. Rinse the membrane with high purity water and equilibrate in transfer buffer for 5 minutes prior to use. Do not allow the membrane to dry out during processing. (If membrane does become partially dry, allow to dry fully, re-wet with 20% alcohol, exchange to buffer, and proceed. This extra alcohol step does not usually interfere with detection.)</p>	<p>Polyvinyliden Fluoride Membrane. Activation: Wet PVDF membrane in methanol for 15 seconds and then transfer to a container filled with ddH<sub>2</sub>O for 5 min. Soak membrane in 1x transfer buffer for 10 min while preparing transfer sandwich. The binding of protein to PVDF is much more sensitive to SDS levels. Too much SDS can inhibit the protein's ability to bind to the PVDF and can, in fact, help proteins already bound to the membrane to slip off. SDS levels should never exceed 0.05% for PVDF.</p>
------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

3. Equilibrate also the gel in transfer buffer (15-20 min).

4. Place 6 absorbent pads on the anode (+) plate.

5. Carefully place the membrane on the saturated pads. Roll a clean glass pipette slowly and gently over the membrane in one direction to eliminate air bubbles that may exist between the pads and the membrane.

6. Place the gel on top of the membrane, rolling a glass pipette slowly and gently over the gel in one direction to eliminate air bubbles that may exist between the gel and membrane.

7. Place 6 absorbent pads on top of the gel, then place the cathode side (-) of the apparatus on top of the stack.

8. Insert the stack in the tank and add transfer buffer per the manufacturer's instructions.

Connect the tank to the power supply and start the transfer. Follow the manufacturer's recommendations for current\*. Transfers are generally complete in 15-90 minutes.

For nitrocellulose membrane: after the transfer is complete, allow the membrane to air dry at room temperature. This helps the proteins to bind more strongly to the membrane, preventing loss during subsequent washes and detection steps.

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### 3. Wet western blot transfer

This protocol describes the procedure for transferring protein from polyacrylamide gel onto membrane using a wet transfer method with a cassette system.


1.

<p>Nitrocellulose membranes should be handled carefully and with gloves. Nitrocellulose membranes can be brittle and may crack or tear, especially when dry. <u>Cut to the size of the SDS-PAGE gel.</u> Slowly lower the membrane into 20% (v/v) ethanol and agitate briefly. The membrane will become translucent as it wets. Rinse the membrane with high purity water and equilibrate in transfer buffer for 5 minutes prior to use. Do not allow the membrane to dry out during processing. (If membrane does become partially dry, allow to dry fully, re-wet with 20% alcohol, exchange to buffer, and proceed. This extra alcohol step does not usually interfere with detection.)</p>	<p>Polyvinyliden Fluoride Membrane. Activation: Wet PVDF membrane in methanol for 15 seconds and then transfer to a container filled with ddH<sub>2</sub>O for 5 min. Soak membrane in 1x transfer buffer for 10 min while preparing transfer sandwich. The binding of protein to PVDF is much more sensitive to SDS levels. Too much SDS can inhibit the protein's ability to bind to the PVDF and can, in fact, help proteins already bound to the membrane to slip off. SDS levels should never exceed 0.05% for PVDF.</p>
------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

2. Equilibrate the gel in transfer buffer (15-20 min).

3. Fill the Blotter tank to about 50% of the fill volume with transfer buffer (0-4°C). (Connect with the cooling system or place an ice block (-20°C) in the dedicated location. Place a magnetic stir bar inside the tank.)

4. Set up the gel membrane sandwich :

	<p>Pour some chilled transfer buffer in the compartment and place the transfer cassette in the compartment of the gel/blot assembly tray. Cathode side (black) is laying horizontal. Place a fiber pad, <u>thoroughly wet</u>. Place a fiber paper piece on top of the fiber pad (it will wet immediately). Place the gel on top of the filter paper. Use roller to remove any air bubbles that may be trapped underneath the gel. Place the membrane on top of the gel taking care not to trap air. The membrane should not be adjusted once it touches the gel. Use roller to roll out bubbles. Place a wet fiber paper piece on top of the membrane. Wet a second fiber pad, <u>completely saturated</u> with transfer buffer. Then place on top of the second filter paper.</p>
------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

5. Lock the cassette.

6. Move the locked cassette into the groove in the blotter tank (red side with red electrode = anode, black side with black electrode = cathode). Add the remaining transfer buffer to fill level. Turn on magnetic agitation.

Connect the tank to the power supply and start the transfer. 100V. Transfers are generally complete in 30-60 minutes.

For nitrocellulose membrane: after the transfer is complete, allow the membrane to air dry at room temperature. This helps the proteins to bind more strongly to the membrane, preventing loss during subsequent washes and detection steps.

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