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200 µL of culture is placed in a 1.5 mL Eppendorf tube. Then celles are lysed by adding 200 µL of 10% sodium dodecyl sulfate, mixing for 5 s with a vortex mixer, and incubating at room temperature for 5 min. Violacein is quantitatively extracted from this cell lysate by adding 900 µL of water-saturated butanol, vortexing for 5 s, and centrifuging at 13 000 g for 5 min in a microfuge. The butanol (upper) phase containing the violacein was collected and its absorbance measured at a wavelength of 585 nm.

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Collection of 0.2 mL of culture. Centrifugation at 10 000 g for 10 min. The supernatant is eliminated and 1 mL absolute ethanol is added to the pellets of the cells to extract the pigment. Vortexing. Centrifugation at 10 000 g for 10 min. The supernatant is possibly diluted in absolute ethanol for violacein (violacein+deoxyviolacein) by absorbance measurement at 575 nm (molar extinction coefficient of pure violacein = 0,05601 ml.µg⁻¹ cm₋₁ in ethanol at 575 nm).

Bibliographie

- Extraction of violacein from *Chromobacterium violaceum* provides a new quantitative bioassay for N-acyl homoserine lactone autoinducers ; J. of Biological Methods (40) 2000 47-55.
- Optimization of culture conditions for violacein production by a new strain of *Duganella* sp. B2 ; Biochemical engineering journal ; 2009, vol. 44, no2-3, pp. 119-124
- Factorial design and response surface optimization of crude violacein for *Chromobacterium violaceum* production ; Biotechnology letter volume 23, décembre 2001, 1963-1969
- <http://web2.uwindsor.ca/courses/biology/fackrell/Microbes/4360.htm>
- Violacein and biofilm production in *Janthinobacterium lividum* ; Journal of Applied Microbiology ; 102 (2007) 992-999
- <http://onlinelibrary.wiley.com/doi/10.1111/j.1365-2672.2006.03155.x/pdf>
- Technical Rule for Biological Agents 466 <http://www.baua.de/en/Topics-from-A-to-Z/Biological-Agents/TRBA/TRBA-466.html> et J Med Assoc Thai. 1992 Mar;75 Suppl 2:6-10 ; Hospital acquired *Janthinobacterium lividum* septicemia in Srinagarind Hospital.
- Claudia Hornung et al ; The *Janthinobacterium* sp. HH01 Genome Encodes a Homologue of the *V. cholerae* CqsA and *L. pneumophila* LqsA Autoinducer Synthases ; Plos one ; Published: February 6, 2013 ; <http://dx.doi.org/10.1371/journal.pone.0055045>
- Natalia Valdes et al ; Draft genome sequence of *Janthinobacterium lividum* strain MTR reveals its mechanism of capnophilic behavior; Standards in Genomic Sciences 2015 10:110 ; <https://www.ncbi.nlm.nih.gov/pubmed/26605004>