## Quorum sensing

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#### ***Chromobacterium violaceum* CV026 as a biosensor:**

*Chromobacterium violaceum* is a Gram-negative bacterium that produces a purple pigment, violacein, under quorum sensing control. The strain used as a biosensor, *Chromobacterium violaceum* CV026, is a mini-*Tn5* mutant strain unable to produce AHL and thus unable to produce the purple pigment. Upon the addition of external AHL (commercial, purified or supernatants that have AHLs), the mutation is complemented and violacein is produced. This is the basis of the direct assay with *Chromobacterium violaceum* CV026, which detects AHLs with acyl-side chains with 4 to 8 carbon atoms. A positive result is the production of a purple halo on a white background.

The presence of antagonistic AHLs can inhibit the complementation of the mutation. Thus, if the biosensor strain is put in the presence of the natural inducer (commercial AHL) and an antagonistic AHL is present as well, the mutation will not be complemented. This is the basis of the reverse assay with *Chromobacterium violaceum* CV026, which detects AHLs with acyl side chains longer than 8 carbon atoms. A positive result is the production of a white halo on a purple background.

(From Dr. Graciela Brelles-Mariño Biological Sciences Department, California State Polytechnic University, 2009.)

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
#### Supernatants preparation

- Pour a few ml of each bacterial suspension to be tested into an tube.
- Centrifuge the tubes at 4000 g (or more) for 5 minutes.
- Remove each supernatant, filter through a sterile 0.2  $\mu\text{m}$  pore size filter and place it into another labeled sterile eppendorf tube.

#### Direct assay

- Seed 3 mL of molten semi-solid medium (LB agar 0,3%(w/v)) at 45 °C (surfusion) with 30  $\mu\text{L}$  of an overnight LB culture of *C. violaceum* CVO26 and pour immediately over the surface of a LB agar plate.
- Divide the plate in 4. Each area will be used to test one bacterial supernatant.
- Let the plate dry.
- Gently pour 10  $\mu\text{L}$  of the supernatant to be tested. Note: Don't break the agar, just pour the drop without touching the surface (or 20  $\mu\text{L}$ , or 40  $\mu\text{L}$ , but it needs a well in the medium)
- In the center of the plate, pour a 5  $\mu\text{L}$  of the commercial AHL C6-HSL (5  $\mu\text{M}$ , N-(3-hexanoyl)-DL-homosérine lactone) as a positive control.
  - Let the plate dry.
  - Incubate the plates at 28 °C, 24 hrs.

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#### Preparation of the reverse assay plates:

- Label a Petri dish with LB agar. Divide the plate in 4. Each area will be used to test one bacterial supernatant.
- Place 5  $\mu\text{L}$  of each bacterial supernatant in each area of the Petri dish. (The *Pseudomonas aeruginosa* supernatant is supposed to be a positive control)
- Let the plate dry.
- While the plate is drying out, prepare the tubes as follows:
- Pour 3 mL of the semi solid molten medium. If you are not going to process the tubes immediately, place them into the water bath. Add 100  $\mu\text{L}$  of a moderately turbid bacterial suspension of CV026 in 3 mL of semi-solid molten LB medium (45°C, surfusion). Vortex the mixture. Add 2  $\mu\text{L}$  of the commercial HHL 5  $\mu\text{M}$  and vortex the mixture. Precaution: the AHL is dissolved in acetonitrile which is a very volatile solvent. Add the commercial AHL into the mixture being sure that the tip touches the mixture.
- Pour this mixture onto the surface of the agar plate prepared as indicated above.
- Let it dry.
- Incubate the plates at 28 °C, 24 hrs.

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**Table 1.** Structures and relative activities of AHL and AHT compounds in the *C. violaceum* CV026 agar plate bioassays

R	X	Chemical name	Abbreviation	Amount applied (nmol)*	
				Activation	Inhibition
CH <sub>3</sub> CH <sub>2</sub> CH <sub>2</sub>	O	N-Butanoyl-L-homoserine lactone	BHL	0-9	NA
CH <sub>3</sub> CH <sub>2</sub> CH <sub>2</sub>	S	N-Butanoyl-L-homocysteine thiolactone	BHT	1-7	NA
CH <sub>3</sub> COCH <sub>2</sub>	O	N-(3-Oxobutanoyl)-L-homoserine lactone	OBHL	3-5	NA
PhCH <sub>2</sub> COCH <sub>2</sub>	O	N-Benzoylacetyl-L-homoserine lactone	BAHL	1-53	NA
CH <sub>3</sub> (CH <sub>2</sub> ) <sub>4</sub>	O	N-Hexanoyl-L-homoserine lactone	HHL	0-03	NA
CH <sub>3</sub> (CH <sub>2</sub> ) <sub>4</sub>	S	N-Hexanoyl-L-homocysteine thiolactone	HHT	0-05	NA
CH <sub>3</sub> (CH <sub>2</sub> ) <sub>2</sub> COCH <sub>2</sub>	O	N-(3-Oxohexanoyl)-L-homoserine lactone	OHHL	0-19	NA
CH <sub>3</sub> (CH <sub>2</sub> ) <sub>2</sub> COCH <sub>2</sub>	O	N-(3-Oxohexanoyl)-D-homoserine lactone	(D)OHHL	12-0	NA
CH <sub>3</sub> (CH <sub>2</sub> ) <sub>2</sub> COCH <sub>2</sub>	S	N-(3-Oxohexanoyl)-L-homocysteine thiolactone	OHHT	0-69	NA
CH <sub>3</sub> (CH <sub>2</sub> ) <sub>6</sub>	O	N-Octanoyl-L-homoserine lactone	OHL	0-22	NA
CH <sub>3</sub> (CH <sub>2</sub> ) <sub>4</sub> COCH <sub>2</sub>	O	N-(3-Oxo-octanoyl)-L-homoserine lactone	OOHL	0-33	NA
CH <sub>3</sub> (CH <sub>2</sub> ) <sub>8</sub>	O	N-Decanoyl-L-homoserine lactone	DHL	NA	0-49
CH <sub>3</sub> (CH <sub>2</sub> ) <sub>6</sub> COCH <sub>2</sub>	O	N-(3-Oxodecanoyl)-L-homoserine lactone	ODHL	NA	4-6
CH <sub>3</sub> (CH <sub>2</sub> ) <sub>10</sub>	O	N-Dodecanoyl-L-homoserine lactone	dDHL	NA	4-4
CH <sub>3</sub> (CH <sub>2</sub> ) <sub>8</sub> COCH <sub>2</sub>	O	N-(3-Oxododecanoyl)-L-homoserine lactone	OdDHL	NA	4-2
CH <sub>3</sub> (CH <sub>2</sub> ) <sub>12</sub>	O	N-Tetradecanoyl-L-homoserine lactone	tDHL	NA	4-0
CH <sub>3</sub> (CH <sub>2</sub> ) <sub>10</sub> COCH <sub>2</sub>	O	N-(3-Oxotetradecanoyl)-L-homoserine lactone	OtDHL	NA	3-8

\*Amount of a given compound added to a well cut in the agar to either induce violacein production (activation) or inhibit HHL-mediated induction of violacein production (inhibition) in CV026. NA, No activity i.e. the compound is unable to either induce or to antagonize HHL-mediated induction of violacein in CV026.

## Bibliographie

- McClean KH, Winson MK, Fish L, Taylor A, Chhabra SR, Camara M, Daykin M, Lamb JH, Swift S, Bycroft BW, Stewart GS, Williams P ; Quorum sensing and *Chromobacterium violaceum* : exploitation of violacein production and inhibition for the detection of N-acylhomoserine lactones ; *Microbiology*. 1997 Dec;143 ( Pt 12):3703-11
- Agrégation de biochimie - génie biologique, épreuve pratique de microbiologie, session 2005
- <http://www.csupomona.edu/~gbrelles/Mic428L/Lab%208.doc>
- Swift, Lynch, Fish, Kirke, Tomás, Stewart and Williams ; Quorum Sensing-Dependent Regulation and Blockade of Exoprotease Production in *Aeromonas hydrophila* ; *Infection and Immunity*, October 1999, p. 5192-5199, Vol. 67, No. 10
- Morohoshi, Ebata, Nakazawa, Kato and Ikeda ; N-acyl Homoserine Lactone-Producing or -Degrading Bacteria Isolated from the Intestinal Microbial Flora of Ayu Fish (*Plecoglossus altivelis*) ; *Microbes and Environments* Vol. 20 (2005) , No. 4 pp.264-268
- Atkinson, Chang, Sockett, Cámara and Williams ; Quorum Sensing in *Yersinia enterocolitica* Controls Swimming and Swarming Motility ; *Journal of Bacteriology*, February 2006, p. 1451-1461, Vol. 188, No. 4
- Jangid, Kong, Patole and Shouche , *uxrI* homologs are universally present in the genus *Aeromonas*, *BMC Microbiol.* 2007; 7: 93  
<http://www.ncbi.nlm.nih.gov/pmc/articles/PMC1367215/>